

Subcritical Physics Testing with the RhoPRO[®] Reactivity Computing System

Background

Prior to beginning a new operating cycle after refueling, Low Power Physics Testing (LPPT) is the first series of tests to demonstrate consistency between the reconstructed core and its nuclear design model. The core design is validated by demonstrating measured-to-predicted agreement of key characteristics such as reactor shutdown capability, core reactivity balance and control, as well as by providing an initial core power distribution check.

For over 40 years, Westinghouse has continuously developed new products and services to save critical path time during reload plant startups, all the while maintaining or enhancing safety and performance. Our Subcritical Physics Testing (SPT) services, which utilize our RhoPRO[®] reactivity computing system, is the latest in that line of innovative startup test methods.

Description

SPT with RhoPRO is a revolutionary program for core design validation that achieves the same objectives as traditional LPPT through a unique and simplified approach. The process consists of performing steady-state comparisons of inverse count rate ratio (ICRR) over a variety of subcritical core configurations and conditions.

SPT with RhoPRO data collection occurs during normal planned subcritical MODE 3 operations (e.g., in parallel with Rod Drop Time Testing) and again coupled with a rod withdrawal approach to criticality in MODE 2. Data collection is seamlessly integrated into pre-existing activities, significantly reducing critical-path testing time. Testing at critical conditions is therefore reduced to only an isothermal temperature coefficient (ITC) measurement, allowing the plant to be pre-configured for power ascension to MODE 1.

SPT is deployed using Westinghouse's state-of-the-art RhoPRO reactivity computer to connect to isolated nuclear instrumentation and temperature signals. Service offerings include use of the RhoPRO system and 24-hour coverage by qualified Westinghouse test engineers.




RhoPRO Reactivity Computer System

Benefits of SPT over LPPT

- SPT is integrated into existing plant startup activities, eliminating the critical path time associated with LPPT. Critical path savings of four hours (or more) can be realized, compared to traditional ANSI/ANS-19.6.1-endorsed methods
- Elimination of other associated LPPT activities (e.g., rod speed adjustments, nuclear instrumentation channel removal and restoration, nuclear heat determination)
- Critical path impact is limited to the ITC measurement, which typically takes less than one hour
- Implementation is possible for any PWR plant
- Human performance and reactivity management benefits include:
 - Data collection conducted at steady-state subcritical conditions
 - No time-sensitive operator actions are required
 - SPT enables “pull-to-critical” startup rather than an infrequently performed dilution to criticality
 - All ex-core detector channels and their associated protective functions remain fully enabled and in service at all times

SPT with RhoPRO[®] Reactivity Computing System

A comparison of the SPT with RhoPRO and typical LPPT evolutions is depicted in the following illustration.



Traditional LPPT is a direct barrier to MODE 1 entry. SPT with RhoPRO simplifies and streamlines post-refueling startup activities between MODE 3 and MODE 1. As a result, core design verification is almost entirely removed from outage critical path, putting the customer back on the grid faster than ever before.

Experience

The Spatially Corrected Inverse Count Rate (SCICR) methodology was originally approved by the U.S. Nuclear Regulatory Commission (NRC) in 2005, which paved the way for release of the Subcritical Rod Worth Measurement (SRWM[™]) application.

The new SPT with RhoPRO application maintains all fundamental benefits of the original SRWM product, while incorporating new features and changes in response to key lessons learned from 50 total SRWM applications at over 15 U.S. units.

- The new generation SPT / SCICR methodology eradicates masking effects by eliminating interaction between measurement and prediction. In doing so, SPT with RhoPRO offers a “pure” measured-to-predicted comparison.

- The SPT with RhoPRO program and results evaluation package enhance overall detectability of core anomalies.
- The improved built-in process controls of SPT optimized the balance between margin of error and sample time.
- The data collection requirements have been significantly reduced as compared to those of the SRWM application. By reducing the necessary data sampling to 1-2 minutes per state point, which is generally within the normal hold time following a reactivity maneuver, a “net zero” critical path impact can be achieved.

The licensing basis of the SCICR methodology, topical report WCAP-16260, Rev. 2, is fully approved by the NRC as of May 2019.

SPT with RhoPRO was successfully utilized in November of 2019 at a U.S. plant.

RhoPRO[®] and SRWM[™] are trademarks or registered trademarks of Westinghouse Electric Company LLC, its affiliates and/or its subsidiaries in the United States of America and may be registered in other countries throughout the world. All rights reserved. Unauthorized use is strictly prohibited. Other names may be trademarks of their respective owners.